




1. Delete Action should be hide for users group and individual users for Admin users – delete option should be available for Controlytics admin.(developer)



Users Group



Individual Users



In user details section remove delete user option

2. JWT Token settings should be available for Controlytics Admin

GRANULES

Security > General

Home
Users
Audit Trail
Settings
Notification center

General policy

Maximum number of failed login attempts, before account is locked
5

In case user account lockout, send notification to email

Password policy

Minimum password length*	Maximum password length
8	72
Minimum number of uppercase letters	Minimum number of lowercase letters
1	1
Minimum number of digits	Minimum number of special characters
0	1
Password expiration period in days	Password reuse frequency in days
90	272

☐ Allow whitespace ☐ Force to reset password if not valid ⓘ

Undo Save

Auto logout settings ⓘ

Auto logout duration (sec)*
901


Undo Save

This auto logout visible for admin. For controlytics admin(developer) other timeout settings should be visible, as there is some restriction like auto logout cannot be less than refresh token, we are facing issue when configuring this time and there is no option for us to change the parameter as it is not visible

3. Admin self-password change issue


When admin changed his password, and login with a new password. Once login success, it asking password confirmation popup sometimes. And also, sometimes password confirmation pops up showing after so much idle time.

We refreshed the webpage also after logout and tried to login, still the same issue came, if we close the popup without entering the password its closing and working as usual.



After admin changed his password show password changed succesfully. popup (Do if possible)


- 4. Attributes values were not proper for Controllytics Admin and delete option should be available for Attributes for controllytics admin.**



For other user's attribute's view was proper only.

5. For Controlytics Admin – when clicking on Audit trail its redirecting to Home page.[DO if time permits]

6. Remove pencil icon in Dashboard list for Supervisor and Operator and Maintenance, Show only Dashboard Names.
Currently Pencil icon is visible for Operator and Supervisor, not visible for Maintenance, It should not be visible for Supervisor,Operator and Maintenance,



7. In Dashboard Time Period Selection, After Selecting Time Period range and updating, its showing YYYY/MM/DD format, its should be DD/MM/YYYY

Block	Area	Location	Instrument ID	Gateway	Temperature	Humidity
CWH	CVC-3	Bottle Packing-IV Secondary Packing	CWH/RDU/076-00	IOTGW03_158	25 °C	
CWH	Module-F	Blender Charging Area-02	CWH/RDU/043-00	IOTGW01_109	18 °C	Error
CWH	Stability Area	Chambers Area	CWH/RDU/025-00	IOTGW01_161	25.2 °C	69.7 %RH

8. Hide Check box selection in Device shared Attributes for Maintenance User

Instrument ID	Instrument Name
☐	Testing Device
☐	CWH/RDU/070-00
☐	CWH/RDU/068-00
☐	CWH/RDU/069-00
☐	CWH/RDU/067-00
☐	CWH/RDU/066-00
☐	CWH/RDU/065-00
☐	CWH/RDU/064-00
☐	CWH/RDU/063-00
☐	CWH/RDU/062-00

Shared attributes	
Key ↑	Value
Area	CWH 1st Floor
Block	CWH
Enable_Humid_NLT	true
Enable_Humid_NMT	true
Enable_Temp_NLT	true
Enable_Temp_NMT	true
High Humidity Limit	80
High Humidity Range	100
High Temperature Limit	27
High Temperature Range	50

9. Edit Option Functionality was not working, so hide the Edit option[pencil] as we dont need it.

 GRANULES

Home

Devices

Functional Locations

Audit Trail

Functional Locations

Functional Location Type

All

☐

Functional Location

☐

CWH Block

☐

Module E

☒

Module F

☐

Microbiology

☐

CVC-IV & III

☐

Stability & IT & Sample Storage

☐

OLD QC

☐

NEW QC

Module F

Functional Location details

Details

Relations

Outbound relations

Direction

From

Type ↑	To entity type	To entity name		
Contains	Device	CWH/RDU/038-00		
Contains	Device	CWH/RDU/033-00		
Contains	Device	CWH/RDU/041-00		
Contains	Device	CWH/RDU/035-00		
Contains	Device	CWH/RDU/039-00		
Contains	Device	CWH/RDU/058-00		
Contains	Device	CWH/RDU/043-00		
Contains	Device	CWH/RDU/040-00		
Contains	Device	CWH/RDU/042-00		
Contains	Device	CWH/RDU/057-00		

Items per page: 10 1 - 10 of 15

10. For auto logout duration configuration - password authentication was not asking and also not recording in Audit Trail.

11. In data report, Change title as Temperature and RH Data Report



Unit Name : Gagillapur

TEMP AND RH LOG REPORT

Block : CWH Block	Start Date : 12-10-2025
Area : Module-E	Start Time : 12:10
Location : Excipient Dispensing Area	End Date : 13-10-2025
Inst ID : CWH/RDU/029-00	End Time : 12:10
Temp. Limit : NMT 27 , NLT 14	RH Limit : NMT 80 , NLT 40
Temp. Max : 24.6 , Temp. Min : 22.8	RH Max : 61.0 , RH Min : 53.0
	Interval : 5 minutes

DATE	TIME	TEMPERATURE (°C)	RELATIVE HUMIDITY(%)
12-10-2025	12:10	N/A	N/A
12-10-2025	12:15	23.5	56
12-10-2025	12:20	23.3	56
12-10-2025	12:25	23.2	56

12. Once user password is expired it's not redirecting to Password Reset Page to change his password.

13. Once Admin received forgot password request, Admin reset user password then its not recording audit trail.

14. Hide Add, refresh, search, check box selection options in Admin role.

